

HYSTEROMETER

A DEVICE TO ASSESS THE TOTAL UTERINE DEPTH

The Hysterometer is a device for preliminary hysterometry.

FEATURES

- three sizes: 10, 12, 14
- 15 depth markings from distal tip
- flexible
- smooth rounded tip for easy insertion






ORIGINAL PICTURE

GYNO INFO HYSTEROMETER

Length 10	284 mm
Length 12	290 mm
Length 14	290 mm

Sterile	EO sterilized
---------	---------------

Packaging	75/shipping box
-----------	-----------------

Order number	#4171-10 
	#4171-12 
	#4171-14 



GMP Belgium Hysterometer 06/09



Gynetics Medical Products n.v.
Rembergt Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine